

# Nimco Test Shoe Feedback and Changes:



- Patient **MUST** be barefoot for best evaluation
- Mark areas of concern with sharpie.
- Evaluate with patient fully weightbearing
- Provide pictures from all angles.
- Tape closed test shoe with masking tape
- Video with verbal explanation if possible.

	Left Good Fit	Right Good Fit	Adjust Fit Left (Please provide detailed measurements)	Adjust Fit Right (Please provide detailed measurements)
Length (option to add pedigraph photo with length measurement changes)	<input type="checkbox"/>	<input type="checkbox"/>	Add _____ mm Subtract _____ mm Notes:	Add _____ mm Subtract _____ mm Notes:
Metatarsal/ Forefoot Width (option to add pedigraph photo with length measurement changes)	<input type="checkbox"/>	<input type="checkbox"/>	Medial Add _____ mm Medial Subtract _____ mm  Lateral Add _____ mm Lateral Subtract _____ mm Notes:	Medial Add _____ mm Medial Subtract _____ mm  Lateral Add _____ mm Lateral Subtract _____ mm Notes:
Volume (This is dictated by the plastic closure. Tape closed and mark either space between two sides or overlap.)	<input type="checkbox"/>	<input type="checkbox"/>	<b>Cone/Instep:</b> Overlap _____ mm Space Between _____ mm <b>Midfoot:</b> Overlap _____ mm Space Between _____ mm <b>Toes Box:</b> Overlap _____ mm Space Between _____ mm Notes:	<b>Cone/Instep:</b> Overlap _____ mm Space Between _____ mm <b>Midfoot:</b> Overlap _____ mm Space Between _____ mm <b>Toes Box:</b> Overlap _____ mm Space Between _____ mm
Depth (With closure taped together as close as possible then measure from the closure to your preferred depth.)	<input type="checkbox"/>	<input type="checkbox"/>	<b>Heel:</b> Add _____ mm Subtract _____ mm <b>Cone/Instep:</b> Add _____ mm Subtract _____ mm <b>Midfoot:</b> Add _____ mm Subtract _____ mm <b>Toes Box:</b> Add _____ mm Subtract _____ mm Notes:	<b>Heel:</b> Add _____ mm Subtract _____ mm <b>Cone/Instep:</b> Add _____ mm Subtract _____ mm <b>Midfoot:</b> Add _____ mm Subtract _____ mm <b>Toes Box:</b> Add _____ mm Subtract _____ mm Notes:
Heel Closure (With test fit taped as close as possible)	<input type="checkbox"/>	<input type="checkbox"/>	Overlap _____ mm Space Between _____ mm Notes:	Overlap _____ mm Space Between _____ mm Notes: